

# Lab Book

Raphael Angerbauer

2024-03-01

# Table of contents

<b>Preface</b>	<b>3</b>
<b>1 Competitive inhibition of transcription using STYX</b>	<b>4</b>
1.1 Flourometry . . . . .	4
1.1.1 Attempt 1 . . . . .	4
1.1.2 Attempt 2 . . . . .	5
1.1.3 Attempt 3 . . . . .	7
1.1.4 Attempt 4 . . . . .	9
1.1.5 Final Result . . . . .	11

# Preface

This is a Laboratory Notebook for experiments conducted from March 2022 to June 2023. The results obtained during this time were used to write the thesis “ ” to obtain the title “Dr. med. univ.”.

The notebook is sorted by experiments and not in chronological order. After each individual repetition has been described, a summary of all the data is shown.

# 1 Competitive inhibition of transcription using STYX

['3-T1', 'Figure4', 'Flourometry', 'WB']

## 1.1 Flourometry

['Results.csv', '1', '2', '3', '4']


### 1.1.1 Attempt 1

#### 1.1.1.1 Seeding (28.11.2022)

- Seed 15.000 HeLa cells per well in 96 Well (diluted in 100µL DMEM)
- Layout:
  - 4 Wells reserved for Blank measurement (termed NT)
  - 4 Wells per row per condition: Conditions: X, TP, FY #### Transfection (29.11.2022)
- Change Media beforehand
- Transfect in a 1:2 manner (Plasmid1:Plasmid2)
- 0.1µg total amount of Plasmid per Well
- Conditions:
  - X: P1 = GFP, P2 = STYX
  - TP: P1 = GFP, P2 = TP-Linker-Protein
  - FY: P1 = GFP, P2 = FY-Linker-Protein
- Change Media 4h after Transfection #### Measurement (30.11.2022)
- Change Media to 100µL PBS
- Measure with Spark @ 488nm
- Incubate on ice for 20mins
- Replace PBS with 10µL RIPA(+PI)
- Incubate for 20mins
- Take out 10µL of cell lysate into new 96-Well

- Put in 150mL of Bradford-assay
- Measure @ 660nm ##### Results The mean value of all 5 measurements in each well is calculated. The intensity of the blank measurement is subtracted. Then the resulting value is divided by the measured protein concentration and then further normalised to the mean of the control.

	mean	std	count
group			
X	1.000	0.043	4
TP	0.549	0.097	4
FY	0.585	0.157	4




### 1.1.2 Attempt 2

#### 1.1.2.1 Seeding (23.08.2022)

- Seed 15.000 HeLa cells per well in 96 Well (diluted in 100µL DMEM)

- Layout:
  - 4 Wells reserved for Blank measurement (termed NT)
  - 4 Wells per row per condition: Conditions: X, TP, FY #### Transfection (24.08.2022)
- Change Media beforehand
- Transfect in a 1:2 manner (Plasmid1:Plasmid2)
- 0.1µg total amount of Plasmid per Well
- Conditions:
  - X: P1 = GFP, P2 = STYX
  - TP: P1 = GFP, P2 = TP-Linker-Protein
  - FY: P1 = GFP, P2 = FY-Linker-Protein
- Change Media 4h after Transfection #### Measurement (25.08.2022)
- Change Media to 100µL PBS
- Measure with Spark @ 488nm
- Incubate on ice for 20mins
- Replace PBS with 10µL RIPA(+PI)
- Incubate for 20mins
- Take out 10µL of cell lysate into new 96-Well
- Put in 150mL of Bradford-assay
- Measure @ 660nm #### Results The mean value of all 5 measurements in each well is calculated. The intensity of the blank measurement is subtracted. Then the resulting value is divided by the measured protein concentration and then further normalised to the mean of the control.

	mean	std	count
group			
X	1.000	0.021	3
TP	0.342	0.023	3
FY	0.409	0.076	3



### 1.1.3 Attempt 3

#### 1.1.3.1 Seeding (05.12.2022)

- Seed 15.000 HeLa cells per well in 96 Well (diluted in 100µL DMEM)
- Layout:
  - 4 Wells reserved for Blank measurement (termed NT)
  - 4 Wells per row per condition: Conditions: X, TP, FY

#### 1.1.3.2 Transfection (06.12.2022)

- Change Media beforehand
- Transfect in a 1:2 manner (Plasmid1:Plasmid2)
- 0.1µg total amount of Plasmid per Well
- Conditions:
  - X: P1 = GFP, P2 = STYX

- TP: P1 = GFP, P2 = TP-Linker-Protein
- FY: P1 = GFP, P2 = FY-Linker-Protein
- Change Media 4h after Transfection

### 1.1.3.3 Measurement (07.12.2022)


- Change Media to 100 $\mu$ L PBS
- Measure with Spark @ 488nm
- Incubate on ice for 20mins
- Replace PBS with 10 $\mu$ L RIPA(+PI)
- Incubate for 20mins
- Take out 10 $\mu$ L of cell lysate into new 96-Well
- Put in 150mL of Bradford-assay
- Measure @ 660nm

### 1.1.3.4 Results

The mean value of all 5 measurements in each well is calculated. The intensity of the blank measurement is subtracted. Then the resulting value is divided by the measured protein concentration and then further normalised to the mean of the control.

	mean	std	count
group			
X	1.000	0.175	4
TP	0.350	0.059	4
FY	0.749	0.073	4





#### 1.1.4 Attempt 4

##### 1.1.4.1 Seeding (06.12.2022)

- Seed 15.000 HeLa cells per well in 96 Well (diluted in 100µL DMEM)
- Layout:
  - 4 Wells reserved for Blank measurement (termed NT)
  - 4 Wells per row per condition: Conditions: X, TP, FY

##### 1.1.4.2 Transfection (07.12.2022)

- Change Media beforehand
- Transfect in a 1:2 manner (Plasmid1:Plasmid2)
- 0.1µg total amount of Plasmid per Well
- Conditions:
  - X: P1 = GFP, P2 = STYX

- TP: P1 = GFP, P2 = TP-Linker-Protein
- FY: P1 = GFP, P2 = FY-Linker-Protein
- Change Media 4h after Transfection


#### 1.1.4.3 Measurement (08.12.2022)

- Change Media to 100 $\mu$ L PBS
- Measure with Spark @ 488nm
- Incubate on ice for 20mins
- Replace PBS with 10 $\mu$ L RIPA(+PI)
- Incubate for 20mins
- Take out 10 $\mu$ L of cell lysate into new 96-Well
- Put in 150mL of Bradford-assay
- Measure @ 660nm

#### 1.1.4.4 Results

The mean value of all 5 measurements in each well is calculated. The intensity of the blank measurement is subtracted. Then the resulting value is divided by the measured protein concentration and then further normalised to the mean of the control.


	mean	std	count
group			
X	1.000	0.188	4
TP	0.461	0.082	4
FY	0.548	0.139	4



### 1.1.5 Final Result

	mean	sd	ci_lo	ci_hi	count
group					
X	1.000	0.117	0.936	1.057	15
TP	0.431	0.106	0.377	0.485	15
FY	0.583	0.156	0.501	0.664	15

	mean [CI]	value	star
group			
X	1.0 [0.936, 1.057]	NaN	NaN
TP	0.431 [0.377, 0.485]	0.0	****
FY	0.583 [0.501, 0.664]	0.0	****



2